



www.medsourcerental.com

MEDSource Inc.

Wide Belt Lead Apron



A fully adjustable 6" wide elastic belt transfers the apron weight from the shoulders to the hips making this apron the most comfortable for longer procedures.Â Just shrug your shoulders and adjust the belt to a comfortable position.Â Shoulder pads offer additional comfort.Â Velcro shoulder strap keeps the apron snug across the upper back so there's no shoulder slippage.Â Pocket is included.



Available in the following lead options (from heaviest to lightest weight):

Cost Cruncher: .35mm and .5mm

True Lite: .35mm and .5mm

Star Lite: .35mm and .5mm

Prestige: .35mm and .5mm

Colors: 24 available colors at no extra charge

Designer Prints: 58 available at an extra charge

Sizes: Small to XL; custom sizes also available for extra charge

Embroidery: Available for extra charge

Â

[Vendor Information](#)