



www.medsourcerental.com

MEDSource Inc.

Veterinary Mittens



Designed for use while positioning animals during procedures, lead veterinary mittens give full wrap-around 0.5 mm lead protection for veterinarians and veterinary assistants. These versatile mittens weigh 4 lbs. and have a Velcro wrist strap for a snug fit.

Unlike palm shields, our mitten design gives full wraparound 0.5mm protection while positioning animals during procedures. velcro wrist strap. NO COLOR CHOICE. 4 lb.

[Vendor Information](#)