



www.medsourcerental.com

MEDSource Inc.

TWC24 Triple Chamber Warming Cabinet



Triple Chamber Warming Cabinet with One Adjustable Shelf

Standard Features: Keyed temperature lock out
Push button digital operation Large dual digital displays Displays actual temperature and set point at all times Separate digitally controlled heating chambers Rapid warm time (2 - 6 hours depending on load) Uniform heating Fully insulated Audible and visual over-temp alarm UL/CUL approved FDA registered OSHPD pre-approved

Additional Options: Full view glass door Mobile dolly or cart 18" depth Full recess Pass through chamber Left hand door hinge Keyed door lock Temperature recorder Seismic brackets 220/240 V, 50/60HZ Custom sizes

*Picture shows OPTIONAL dolly (additional fee)

[Vendor Information](#)

Overall Size

30" W x 26.5" D x 78.25" T

14" cubic feet of storage

1" adjustable shelf **Upper Chamber I.D.** 26" W x 23" D x 11.25" T 3.9 cubic feet of storage

Temperature range: 90-110 F. Accuracy +/- 1 degree F. of set point **Intermediate Chamber I.D.** 26" W x 23" D x 10.75" T 3.7 cubic feet of storage
Temperature range: 90-135 F. Accuracy +/- 1 degree F. of set point

Lower Chamber I.D. 26" W x 23" D x 18.5" T 6.4 cubic feet of storage
Temperature range: 90-160 F. Accuracy +/- 1 degree F. of set point
Electrical 120V, 13.5 amp., 50/60HZ .65 kw/hr, 2238 btu/hr UL/CUL listed FDA registered